

HD-350 VIDEO PROCESSOR

HD video outputs	DVI, VGA, SDI
Standard video outputs	S-Video, CVBS
HD video recording	Built-in workstation, timely record and review
Zoom / PBP	Yes / Yes
Dimension (L x W x H / weight )	455mm x 370mm x 124mm / 9.5 kg

HDL-35E LIGHT SOURCE

Main lamp	50W LED
Transparency function	7 seconds maximum brightness flicker
Lifetime indicator	Yes
Hot plugging	Endoscope switch without powering off
Dimension (L x W x H / weight )	460mm x 370mm x 154mm / 11kg

430 SERIES SCOPES

Parameter	Gastroscope		Colonoscope			
	EG-430	EG-430L	EC-430	EC-430T	EC-430L	EC-430L/T
Depth of view (mm)	3-100		3-100			
Field of view	140°		140°			
Distal end diameter (mm)	9.3	9.8	12	12	12.9	12.9
Insertion diameter(mm)	9.3	9.8	12.5	12.5	12.9	12.9
Instrument channel diameter (mm)	2.8	3.2	3.8		4.2	
Angulations	UP 210°/DOWN 90° LEFT 100°/RIGHT 100°		UP 180°/DOWN 180° LEFT 160°/RIGHT 160°			
Working length (mm)	1050		1350	1700	1350	1700
Total length (mm)	1400		1700	2050	1700	2050

HD-350

SEE MORE,  
ACCOMPLISH MORE,  
BENEFIT MORE.



Representante exclusivo:

ENDOMED

Tecnologies S.A.C.

Dirección : Av. Los Álamos 309, Urb. José Olaya  
San Juan de Miraflores


Teléfono : (01) 368 - 7542

E-mail : endomedsac@gmail.com

Web : www.endomed.com.pe



E-HD35020240812




Innovation rewrites the standard. With that in mind, SonaScape Medical presents you the HD-350, a brand new high definition endoscopy system designed specifically for screening and diagnosis with an excellent budget friendliness.



## CMOS

With the CMOS image sensor, HD-350 gives high quality images with low noise levels, which leads to improved clarity and optimized vision in not only routine procedures but also in complex interventions.



## Smart Workflow

To increase efficiency, patient data management is accessible with the press of a button. The smart workflow of the built-in workstation of HD-350 helps doctors to edit reports, record videos and manage patient files with ease.



## Clean and Clear

The waterjet function of the scopes helps to maintain a clear endoscopic view by washing away blood or mucous.

